SPANISH EDITION

**Motec<sup>®</sup>**

Wrist Joint Prosthesis

*Swemac*

# Motec®

## Wrist Joint Prosthesis

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“

*... rotation occurs about a fixed axis located within the head of the capitate, and the location of each axis is not changed by the position of the hand in either plane.*

”

Youm Y, McMurthy RY, Flatt AE, Gillespie TE.

An experimental study of radial-ulnar deviation and flexion-extension.

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
*CFR-PEEK represents an alternative load-bearing material because of its superior mechanical and chemical behaviour without any increased biological activity of the wear particles, compared with a standard load-bearing material.*

Utzschneider S, Becker F, Grupp TM, Sievers B, Paulus A, Gottschalk O, Jansson V.

Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo.

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## Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients

Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.


*J Hand Surg Am. 2017*

### Abstract

Fifty-seven (40 male) patients with end-stage arthritis changes received an uncemented ball-and-socket total wrist arthroplasty (Motec Wrist). Function was evaluated before surgery and at yearly follow-ups. Visual analog scale at rest and activity, quick Disabilities of the Arm, Shoulder, and Hand (QuickDASH), active range of motion (AROM), and grip-strength were recorded. Standardized radiographs were taken to assess osteolysis, loosening, and subsidence.

Fifty-six patients were followed for a mean of 8 years (SD, 2 years). Eight wrists were reoperated with arthrodesis (4) or a new arthroplasty (4) owing to distal component loosening (3), infection (2), pain/fixed malposition (2), or proximal and distal component loosening (1). One radiocarpal dislocation was reduced closed and remained stable. Improved QuickDASH score and visual analog scale pain score both at rest and during activity were found at the last follow-up, as well as increased AROM (97 vs 126) and grip strength (21 kg vs 24 kg). The radiological follow-up demonstrated loosening in 2 wrists. Thirty-five patients were working at surgery (17 manual labor) and 27 (11 manual labor) at follow-up. The 10-year Kaplan-Meier survival of the implants was 86% for revision any cause, 2 additional arthroplasties are loose (but not revised), giving a survival rate of 82% if these are revised prior to 10 years of observation.

An uncemented total wrist arthroplasty can provide long-lasting unrestricted hand function in young and active patients.

Kaplan-Meier survival curve

# Referencias

1. **Reigstad O, Holm-Glad T, Bolstad B, Grimsgaard C, Thorkildsen R, Røkkum M.** Five- to 10-Year Prospective Follow-Up of Wrist Arthroplasty in 56 Nonrheumatoid Patients. *J Hand Surg Am.* 2017
2. **Glad TH, Røkkum M, Thorkildsen R, Reigstad O.** Rearticulation from arthrodesis with Motec prosthesis (Reartikulering av avstivede håndledd med Motec protese). *NOF Høstmøteboken Abstract 210*, 2016.
3. **Youm Y, McMurthy RY, Flatt AE, Gillespie TE.** Kinematics of the wrist. I. An experimental study of radial-ulnar deviation and flexion-extension. *J Bone Joint Surg Am.* 1978 Jun;60(4):423-31.
4. **Utzschneider S, Becker F, Grupp TM, Sievers B, Paulus A, Gottschalk O, Jansson V.** Inflammatory response against different carbon fiber-reinforced PEEK wear particles compared with UHMWPE in vivo. *Acta Biomater.* 2010 Nov;6(11):4296-304.
5. **Reigstad A, Reigstad O, Grimsgaard C, Røkkum M.** New concept for total wrist replacement. *J Plast Surg Hand Surg.* 2011 Jun;45(3):148-56.
6. **Glad TH, Røkkum M, Thorkildsen R, Reigstad O.** Chrome and Cobalt blood levels in a metal-on-metal total wrist arthroplasty. *NOF Høstmøteboken Abstract 212*, 2016.
7. **Lundborg G, Besjakov J, Brånemark PI.** Osseointegrated wrist-joint prostheses: a 15-year follow-up with focus on bony fixation. *Scand J Plast Reconstr Surg Hand Surg.* 2007;41(3):130-7.
8. **Røkkum M.** On Late Complications With Ha Coated Hip Arthroplasties, Department of Biomaterials/Handicap Research, Institute for Surgical Sciences, Faculty of Medicine, University of Göteborg, Göteborg, Sweden and Orthopaedic University Clinic, National Hospital, Oslo, Norway, Göteborg 2001.
9. **Reigstad O, Franke-Stenport V, Johansson CB, Wennerberg A, Røkkum M, Reigstad A.** Improved bone ingrowth and fixation with a thin calcium phosphate coating intended for complete resorption. *J Biomed Mater Res B Appl Biomater.* 2007 Oct;83(1):9-15.
10. **Reigstad O, Johansson C, Stenport V, Wennerberg A, Reigstad A, Røkkum M.** Different patterns of bone fixation with hydroxyapatite and resorbable CaP coatings in the rabbit tibia at 6, 12, and 52 weeks. *J Biomed Mater Res B Appl Biomater.* 2011 Oct;99(1):14-20.
11. **Reigstad O, Lütken T, Grimsgaard C, Bolstad B, Thorkildsen R, Røkkum M.** Promising one- to six-year results with the Motec wrist arthroplasty in patients with post-traumatic osteoarthritis. *J Bone Joint Surg Br.* 2012 Nov;94(11):1540-5.
12. **Reigstad O, Røkkum M.** Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg *J Wrist Surg* 2014;3:211–215
13. **Scholes SC, Unsworth A.** Pitch-based carbon-fibre-reinforced poly (ether-ether-ketone) OPTIMA assessed as a bearing material in a mobile bearing unicondylar knee joint. *Proc Inst Mech Eng H.* 2009 Jan;223(1):13-25.
14. **Scholes SC, Unsworth A.** Wear studies on the likely performance of CFR-PEEK/CoCrMo for use as artificial joint bearing materials. *J Mater Sci Mater Med.* 2009 Jan;20(1):163-70.
15. **Kabir K, Schwiesau J, Burger C, Pflugmacher R, Grupp T, Wirtz DC.** Comparison of Biological Response to UHMWPE and CFR-PEEK Particles in Epidural Space. *Universitätsklinikum Bonn, Department for Orthopaedics and Trauma Surgery, Bonn, Germany, Aesculap AG Research & Development, Tuttlingen, Germany.* 2011.
16. **Grupp TM, Utzschneider S, Schröder C, Schwiesau J, Fritz B, Maas A, Blömer W, Jansson V.** Biotribology of alternative bearing materials for unicompartamental knee arthroplasty. *Acta Biomater.* 2010 Sep;6(9):3601-10. Erratum in: *Acta Biomater.* 2012 Apr;8(4):1659.
17. **Wennerberg A.** On surface roughness and implant incorporation. *Department of Biomaterial/Handicap Research, Göteborg, Sweden.* 1996.

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“

*Conversion of a failed total wrist arthroplasty to arthrodesis can be difficult. A custom made titanium alloy peg was constructed to enable arthrodesis with the original arthroplasty components in situ. Two out of three patients were especially challenging cases with little bone available. Bony union was achieved in 2 to 3 months. The peg simplified a difficult revision situation and gave good, predictable results at follow-up.*

”


Reigstad O, Røkkum M.

Conversion of Total Wrist Arthroplasty to Arthrodesis with a Custom-Made Peg  
J Wrist Surg 2014;3:211–215

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







# Product information








- Needed for CFR-PEEK articulation
- Needed for CoCrMo articulation

## Implants










|                                                    |          |     |                                                                                       |
|----------------------------------------------------|----------|-----|---------------------------------------------------------------------------------------|
| Radius Cup   CoCrMo   Ø15 mm                       | 40-1015S | ●   |    |
| Radius Cup   CoCrMo   Ø18 mm (optional)            | 40-1018S | ●   |    |
| Radius Cup   CFR-PEEK   Ø15 mm                     | 40-1915S | ●   |    |
| Metacarpal Head   Ø15 mm   Short Neck (optional)   | 40-1115S | ● ● |    |
| Metacarpal Head   Ø18 mm   Short Neck (optional)   | 40-1118S | ●   |    |
| Metacarpal Head   Ø15 mm   Medium Neck             | 40-1715S | ● ● |    |
| Metacarpal Head   Ø18 mm   Medium Neck (optional)  | 40-1718S | ●   |   |
| Metacarpal Head   Ø15 mm   Long Neck               | 40-1215S | ● ● |  |
| Metacarpal Head   Ø18 mm   Long Neck (optional)    | 40-1218S | ●   |  |
| Metacarpal Head   Ø15 mm   Extra Long Neck         | 40-1315S | ● ● |  |
| Radius Threaded Implant   length 32 mm             | 40-1332S | ● ● |  |
| Radius Threaded Implant   length 38 mm             | 40-1338S | ● ● |  |
| Radius Threaded Implant   length 44 mm             | 40-1344S | ● ● |  |
| Radius Threaded Implant   length 50 mm             | 40-1350S | ● ● |  |
| Metacarpal Threaded Implant   length 45 mm   Large | 40-1445S | ● ● |  |
| Metacarpal Threaded Implant   length 50 mm   Large | 40-1450S | ● ● |  |
| Metacarpal Threaded Implant   length 55 mm   Large | 40-1455S | ● ● |  |
| Metacarpal Threaded Implant   length 60 mm   Large | 40-1460S | ● ● |  |

|                                                               |          |     |                                                                                     |
|---------------------------------------------------------------|----------|-----|-------------------------------------------------------------------------------------|
| Metacarpal Threaded Implant   length 65 mm   Large (optional) | 40-1465S | ● ● |  |
| Metacarpal Threaded Implant   length 70 mm   Large (optional) | 40-1470S | ● ● |  |
| Metacarpal Threaded Implant   length 45 mm   Small            | 40-1475S | ● ● |  |
| Metacarpal Threaded Implant   length 50 mm   Small            | 40-1480S | ● ● |  |
| Metacarpal Threaded Implant   length 55 mm   Small            | 40-1485S | ● ● |  |
| Metacarpal Threaded Implant   length 60 mm   Small            | 40-1490S | ● ● |  |
| Metacarpal Threaded Implant   length 65 mm   Small (optional) | 40-1495S | ● ● |  |
| Metacarpal Threaded Implant   length 70 mm   Small (optional) | 40-1400S | ● ● |  |

## Trials

|                                                           |         |     |                                                                                       |
|-----------------------------------------------------------|---------|-----|---------------------------------------------------------------------------------------|
| Trial – Radius Cup   Ø15 mm                               | 40-1522 | ●   |  |
| Trial – Radius Cup   Ø18 mm (optional)                    | 40-1521 | ●   |  |
| Trial – Radius Cup   Ø15 mm   For CFR-PEEK Cup            | 40-1541 | ●   |  |
| Trial – Metacarpal Head   Ø15 mm   Short Neck (optional)  | 40-1529 | ● ● |  |
| Trial – Metacarpal Head   Ø18 mm   Short Neck (optional)  | 40-1527 | ●   |  |
| Trial – Metacarpal Head   Ø15 mm   Medium Neck            | 40-1524 | ● ● |  |
| Trial – Metacarpal Head   Ø18 mm   Medium Neck (optional) | 40-1523 | ●   |  |
| Trial – Metacarpal Head   Ø15 mm   Long Neck              | 40-1528 | ● ● |  |
| Trial – Metacarpal Head   Ø18 mm   Long Neck (optional)   | 40-1526 | ●   |  |
| Trial – Metacarpal Head   Ø15 mm   Extra Long Neck        | 40-1602 | ● ● |  |

## Instruments

|                                                        |         |     |                                                                                       |
|--------------------------------------------------------|---------|-----|---------------------------------------------------------------------------------------|
| Hohmann Retractor                                      | 40-1503 | ● ● |    |
| Bits 3,5 mm HEX with Quick-Lock                        | 40-1513 | ● ● |    |
| Impactor                                               | 40-1516 | ● ● |    |
| Guide Wire T-handle                                    | 40-1518 | ● ● |    |
| Cup Remover                                            | 40-1519 | ● ● |    |
| Cannulated Drill for Radius   32-50 mm                 | 40-1546 | ● ● |    |
| Cannulated Drill for Metacarpal III   45-70 mm   Large | 40-1551 | ● ● |    |
| Cannulated Drill for Metacarpal III   45-70 mm   Small | 40-1552 | ● ● |   |
| Guide Wire with sharp tip   Ø2 mm                      | 40-1561 | ● ● |  |
| Guide Wire with round tip   Ø2 mm                      | 40-1563 | ● ● |  |
| Radius Spherical Drill   Ø18 mm (optional)             | 40-1566 | ●   |  |
| Radius Spherical Drill   Ø15 mm                        | 40-1567 | ● ● |  |
| Handle Tri-Lobe with Quick-Lock                        | 45-2585 | ● ● |  |
| Handle Tri-Lobe with Ratchet (optional)                | 40-2593 | ● ● |  |
| Adapter, from AO male to Tri-Lobe female (optional)    | 40-5000 | ● ● |  |
| Hammer                                                 | 52-2211 | ● ● |  |
| Awl                                                    | 23.4997 | ● ● |  |
| Tray and lid                                           | 40-1600 | ● ● |                                                                                       |

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Swemac develops and promotes innovative solutions for fracture treatment and joint replacement. We create outstanding value for our clients and their patients by being a very competent and reliable partner.

# Swemac

Motec Wrist Joint Prosthesis

Fabricante:

**Swemac Innovation AB**

CE 0413

Cobolgatan 1 • SE-583 35 Linköping • Suecia

+46 13 37 40 30 • [info@swemac.com](mailto:info@swemac.com) • [www.swemac.com](http://www.swemac.com)

P125-28-1-ES-20180615

Release date: 2018-06-15

This brochure is a translation of: P125-28-1-20180605