



- Convenient instrumentation and straight forward technique to harvest cancellous autograft bone
- Available as a small, easy to use set compatible with the Gorilla® R3CON Plating System



and that the took oquit

SYSTEM CONTENTS

Bone Graft Harvest System




Bone Graft Harvest System ANCILLARY PRODUCTS

BONE GRAFT HARVEST SYSTEM

PRESERVE™ Bone Wedge System

- Allograft Features:
 - Aseptically processed
 - No gamma irradiation or bleach
 - Density matched
 - Available for Evans, Cotton, MTP, Lapidus, Subtalar, Calc-Cuboid and Hammertoe procedures



V92[™] and V92 FC[™] Cellular Bone Matrix

Allogenic bone matrix containing osteoprogenitor cells which aid in bone remodeling



- 92% Cell Viability Post Thaw
- Final handling properties yield a cohesive, wet sand consistency with V92[™], and a moldable paste consistency with V92 FC[™]
- 150,000 Viable Cells/cc Available for Transplantation Primitive, osteogenic cells only
- Novel Non-DMSO Cryoprotectant
 - No decanting required
- Optimized Bone Scaffold Component
 - Mineralized scaffold is a combination of cancellous and cortical bone
- Safe Tissue Processing
 - All graft portions from a single donor
 - Aseptically processed

CELLULAR SOTS ACTUAL

OTHER PARAGON 28 BIOLOGIC SOLUTIONS





HARV-01 RevA

[™]Trademarks and ® Registered Trademarks of Paragon 28, Inc. © Copyright 2018 Paragon 28, Inc. All rights reserved. Patents: www.paragon28.com/patents

Paragon 28, Inc. 48 Inverness Ct. E., Suite 280 Englewood, CO 80112 USA (855) 786-2828

Paragon 28 Medical Devices Trading Limited 43 Fitzwilliam Square West Dublin 2, D02 K792, Ireland +353(0) 15414756

For additional information on Paragon 28® and its products please visit WWW.paragon28.com

